

Nitinol Tipless Baskets

HALO™ 1.5 Fr Sheath Diameter



- Exceptional irrigation enhances visibility.
- Enhanced kink resistance increases longevity.
- Can share 3.6 Fr channel with a laser fiber.

6-Wire HALO™ 1.5 Fr Sheath Diameter

- Effectively captures both small and large stones.
- Able to extract multiple stone fragments simultaneously.

Vantage™ 2.4 Fr Sheath Diameter


- Maintains basket symmetry.
- Enhances capture of small fragments.
- Minimizes tissue trauma.
- Increased rigidity for larger stone fragments.



Helical Baskets

HERCULES™ Stainless Steel Paired Helical Basket

- Maximum dilation force.
- Available with or without filiform tip.




Flat Wire Baskets

APEX™ Nitinol Flat Wire Basket

- Maximum basket durability.

Summit™ Stainless Steel Flat Wire Basket

- Maximum dilation force.
- Always includes filiform tip.



Handles

- Allows independent basket rotation.
- Ultra lightweight and highly tactile.

Summit does not have handle with rotation feature.



Nitinol Tipless Baskets

| ORDER NO. | SHEATH DIA. (Fr) | SHEATH LENGTH (cm) | SHEATH MATERIAL | BASKET DIA. (mm) | NO. OF WIRES | TIP | COMMENTS |
|---------------------|-------------------------------|--------------------|-----------------|------------------|--------------|---------|---------------------------|
| HALO™ | Nitinol Tipless Basket | | | | | | Optimal Irrigation |
| 554120 | 1.5 | 120 | Polyimide | 12 | 4 | Tipless | |
| 6-Wire HALO™ | Nitinol Tipless Basket | | | | | | Optimal Irrigation |
| 556160 | 1.5 | 120 | Polyimide | 16 | 6 | Tipless | |
| Vantage™ | Nitinol Tipless Basket | | | | | | |
| 550310 | 2.4 | 90 | Polyimide | 16 | 4 | Tipless | |
| 550320 | 2.4 | 120 | Polyimide | 16 | 4 | Tipless | |

Helical Baskets

| ORDER NO. | SHEATH DIA. (Fr) | SHEATH LENGTH (cm) | SHEATH MATERIAL | BASKET DIA. (mm) | NO. OF WIRES | TIP | COMMENTS |
|------------------|--|--------------------|-----------------|------------------|-------------------|--------------|-------------------------------|
| HERCULES™ | Stainless Steel Paired Helical Basket | | | | | | Maximum Dilation Force |
| 530520 | 3.0 | 120 | Polyimide | 11 | 2 Singles, 2 Pair | "O" Tip | |
| 530525 | 3.0 | 120 | Polyimide | 11 | 2 Singles, 2 Pair | 5cm filiform | |

Flat Wire Baskets

| ORDER NO. | SHEATH DIA. (Fr) | SHEATH LENGTH (cm) | SHEATH MATERIAL | BASKET DIA. (mm) | NO. OF WIRES | TIP | COMMENTS |
|----------------|---|--------------------|-----------------|------------------|--------------|--------------|-------------------------------|
| APEX™ | Nitinol Flat Wire Basket | | | | | | Kink Resistance |
| 544310 | 2.4 | 90 | Polyimide | 15 | 4 | "O" Tip | |
| 544320 | 2.4 | 120 | Polyimide | 15 | 4 | "O" Tip | |
| Summit™ | Stainless Steel Flat Wire Basket | | | | | | Maximum Dilation Force |
| 342515 | 3.0 | 90 | Polyimide | 20 | 4 | 5cm filiform | |

Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician.
U.S. Patents: 6,190,394; 8,021,372; 8,523,879 and D457955. Patents pending.

Sacred Heart Medical, Inc.

Distributed By:

6018 Blue Circle Drive
Minnetonka, MN 55343 USA

Tel: (952) 345-0280
Fax: (952) 942-7590
contact@sacredheartmedical.com
www.sacredheartmedical.com

©2016 Sacred Heart Medical, Inc.
All rights reserved. Rev. 03-2016